INSTITUTE FOR MEDICAL RESEARCH KUALA LUMPUR

Performance of Cellex qSARS-CoV-2 IgG/IgM Cassette Rapid Test

Intended use

For the qualitative detection and differentiation of IgM and IgG antibodies to SARS-CoV-2 in serum, plasma (EDTA, citrate) or venipuncture whole blood specimens from patients suspected of COVID-19 infection.

Manufacturar

Cellex, Inc. (United States of America)

Test Principle

A lateral flow immunoassay

Test Kit

The evaluation was carried out using this kit with the lot number of ZCFX20200609WI5513C-3, and the expiry date was Dec 05, 2022.

Instrument Used

None

Reagent and Sample Preparation, Result Interpretation

Kindly refer to product package insert in the attachment

Sample Used

Known Positive Sample= 50 Known Negative Sample= 50 Total samples used for analysis = 100

Performance Analysis

Test		(it Assay	Interpretation
	Positive	Negative	
SARS-CoV2 Positive (COVID-19)	49	1	Sensitivity= 98%
SARS-CoV2 Negative (COVID-19)	0	50	Specificity = 100%
		SARS-CoV2 Positive (COVID-19) Positive 49	SARS-CoV2 Positive (COVID-19) 49 1

Comments

None

Test peri	formed by	Reviewed by	Approved by
Oign			Al alvoro
Pn. Nur Izwani Shaiful Bahrin Evaluator 1 / Officer in charge	Dr. Rozainenee Mohd Zain Evaluator 2 / Officer in charge	Dr. Ravindran A/L Thayan Head of Virology Unit, IDRC, IMR	Dr. Hj. Tahir bin Aris Director of IMR

Disclaimer:

The content of this report does not imply that this product is endorsed or recommended by the Institute for Medical Research (IMR). This report may not be reviewed, abstracted, quoted, reproduced or translated, in part or in whole, without the prior written permission of IMR. IMR does not allow the use of this report or reference in any manner in the labeling or advertising of this product or kit package. However, the permission to import, export, or place in the market is subjected to Medical Device Authority, Ministry of Health Malaysia's approval.